Instructions for Use – EXEL Disposable Syringes

REF	Model
26048	EXEL Disposable Syringe, 1mL Luer Slip Tip w/ LDS Gasket
26049	EXEL Disposable Syringe, 1mL Luer Lock Tip w/ LDS Gasket
26050	EXEL Disposable Syringe, 1mL Luer Lock Tip
26200	EXEL Disposable Syringe, 3mL Luer Lock Tip
26201	EXEL Disposable Syringe, 3mL Luer Slip Tip
26230	EXEL Disposable Syringe, 5mL Luer Lock Tip
26231	EXEL Disposable Syringe, 5mL Luer Slip Tip
26265	EXEL Disposable Syringe, 10mL Luer Lock Tip
26266	EXEL Disposable Syringe, 10mL Luer Slip Tip
26280	EXEL Disposable Syringe, 20mL Luer Lock Tip
26281	EXEL Disposable Syringe, 20mL Luer Slip Tip
26290	EXEL Disposable Syringe, 30mL Luer Lock Tip
26291	EXEL Disposable Syringe, 30mL Luer Slip Tip, Eccentric
26300	EXEL Disposable Syringe, 50mL Luer Lock Tip
26301	EXEL Disposable Syringe, 50mL Luer Slip Tip, Eccentric
26305	EXEL Disposable Syringe, 50mL Luer Slip Tip

1 Device Description

- 1.1 EXEL Disposable Syringes are piston syringes of a three-piece design consisting of a calibrated hollow barrel, a plunger, and a gasket. At the tip of the barrel there is a male connector (nozzle) for fitting the female connector (hub) of a hypodermic single lumen needle or other compatible device. The types of nozzle are as follows:

- 1.1.1 **Luer Lock:** Secures a needle onto the syringe by locking the needle hub onto the nozzle. The luer-lock makes a stable connection between the syringe and needle.
- 1.1.2 **Luer Slip:** Secures the needle only with a compression fitting of the needle hub to the syringe nozzle.
- 1.1.3 **Luer Slip, Eccentric:** A luer slip nozzle that is off-center (almost flush with the side of the barrel) for accessing injection points that are difficult to reach with a centered nozzle.

The syringes are packaged with a protective cap over the nozzle which is removed before use.

- 1.2 The volume of solution inside a syringe is indicated by the graduation lines on the barrel. Graduation lines may be indicated in milliliters or fractions of a milliliter, depending on the model.

- 1.2.1 EXEL syringes are labeled based on their nominal capacity (e.g., 3mL, 5mL, 10mL, 50mL), which corresponds with the industry standard specification for each size. Some syringes include extended graduation marks beyond the nominal capacity:

- 5mL syringes are graduated up to 6mL of total graduations
- 10mL syringes are graduated up to 12mL of total graduations
- 20mL syringes are graduated up to 25mL of total graduations
- 30mL syringes are graduated up to 35mL of total graduations
- 50mL syringes are graduated up to 60mL of total graduations

The extended portion is differentiated with smaller graduation marks and may be used to provide measurement beyond the nominal capacity as needed. Measurement accuracy is guaranteed and validated throughout the entire graduated range, including the extended graduations.

1.3 Materials:

- 1.3.1 Barrel/Nozzle/Plunger/Cap: PP (Ethylene-Polypropylene copolymer)
- 1.3.2 Gasket: Polyisoprene synthetic rubber
- 1.3.3 Lubricant: Silicone oil

2 Indications for Use

- 2.1 EXEL Disposable Syringes are intended to be used:
- 2.1.1 to inject fluids into, or withdraw fluids from, the body.
 - 2.1.2 to administer and/or withdraw (aspirate) fluids/drugs.

3 Intended Users / Patient Population

- 3.1 Intended Users: Qualified healthcare professionals in a professional healthcare environment.
- 3.2 Intended Patient Population: The general population, including both adult and pediatric patients of all age groups, as determined appropriate by the qualified healthcare professional using the device.

4 Clinical Benefit

- 4.1 Safe, easy to use.
- 4.2 Rapid delivery of fluids.
- 4.3 Single-use disposables prevent cross-contamination between patients.

5 Warnings / Precautions

- 5.1 Single use only – Use once and discard immediately in accordance with local safety regulations.
- 5.2 Not intended for drug storage.
- 5.3 Do not reuse. Reuse of these devices can result in contamination or transmit disease.
- 5.4 Do not resterilize or reprocess.
- 5.5 These devices are not intended for intraocular use. Use of these devices for intraocular administration may cause adverse effects due to the presence of silicone oil.
- 5.6 To prevent contamination, dispose of the device immediately following use in a designated biohazard container in accordance with OSHA and CDC guidelines and/or other applicable local safety regulations.
- 5.7 U.S. Federal Law restricts this device to sale by or on the order of a physician (or other properly licensed practitioner).

6 Contraindications

- 6.1 There are no recognized contraindications for the use of these devices.

7 Sterility

- 7.1 EXEL Disposable Syringes are sterilized by EO gas sterilization.
- 7.2 Single use only. Do not reuse.
- 7.3 Do not resterilize or reprocess.
- 7.4 Inspect each package prior to use. Do not use the device if any seal or cavity is damaged or breached or if the expiration date has been exceeded. Once the sterile packaging has been opened the device must be immediately used or discarded.

8 Operating Instructions

- 8.1 Open the sterile packaging by separating at the upper package edge and peeling apart to the point below the nozzle.
- 8.2 Prepare the insertion site with appropriate antiseptic as applicable.
- 8.3 Remove the protective cap from the syringe nozzle and affix to needle hub (or other compatible female luer) using aseptic technique.
 - 8.3.1 For Luer Lock Tip: Push the needle hub or female luer onto the tip and twist clockwise until it is fully seated and secure.
 - 8.3.2 For Luer Slip Tip or Luer Slip Tip, Eccentric: Push the needle hub or female luer straight onto the tip until it is securely in place.

- 8.4 Remove the needle cap (if applicable) and perform the procedure according to established technique. Pull back on the plunger to fill the syringe with fluid and push down on the plunger to express the fluid and empty the syringe.

- 8.5 If injecting fluid into the human body, check for air bubbles prior to injection to reduce the risk of air embolisms:

- 8.5.1 Hold the syringe with the tip pointed upwards.

- 8.5.2 Tap on the syringe. If air bubbles appear, tap until they rise to the top, then push the plunger until a small amount of liquid is expelled.

- 8.6 After the procedure is complete, dispose of the syringe in an appropriate container in accordance with OSHA and CDC guidelines and/or other applicable local safety regulations regarding medical waste.

9 Risks / Side Effects

- 9.1 Use of these devices may expose the user or patient to the following risks and/or side effects:

- 9.1.1 Using the syringes to store drugs or medication can lead cause contamination. Do not use for drug storage.

- 9.1.2 Reuse of these devices may cause transmission of disease or bloodborne pathogens due to contamination. Do not reuse.

- 9.1.3 Use of these devices for intraocular applications can expose patients to risks due to the presence of silicone oil, including persistent vitreous floaters, inflammation, and increased ocular pressure.

- 9.1.4 The presence of air bubbles can expose patients to risk of air embolisms. Ensure that all air is expelled from the syringe before injection.

10 Storage

- 10.1 Keep away from direct sunlight.
- 10.2 Keep dry.

NOTICE

© 2026 Exelint International, Co. All rights reserved. EXEL and the EXEL logo are trademarks of EXELINT International, Co.

Any serious incident that occurs must be reported to the manufacturer and competent authority according to local regulations. EXELINT International, Co. disclaims any responsibility for possible consequences resulting from reuse, resterilization or improper use.

EXELINT International, Co.
2500 Santa Fe Avenue,
Redondo Beach, CA 90278 U.S.A.
Tel: (310) 649-0707 Fax: (310) 649-1178
e-mail: info@exelmed.com

FOR U.S.A. INFO CALL: 800-940-3935
www.exelint.com
e-mail: info@exelint.com

EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

EMERGO AUSTRALIA
Level 20 Tower II Darling Park
201 Sussex Street
Sydney, NSW 2000 Australia

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

Graphic Symbol Legend for Medical Device Labeling

	Medical Device
	Catalog number
	Federal law (USA) restricts this device to sale by or on the order of a physician.
	Consult instructions for use
	Caution
	Manufacturer's Lot Number
	Date of Manufacture
	Unique Device Identifier
	Use by date
	Sterilized by Ethylene Oxide. Sterility guaranteed in unopened and undamaged packaging.
	Sterile packaging (single barrier system)
	Protective packaging outside of sterile barrier
	Non-Pyrogenic
	Single use only. Do not reuse.
	Do not resterilize.
	Do not use if package is damaged
	Keep away from sunlight
	Keep Dry
	This device is not manufactured with Natural Rubber Latex
	Manufacturer
	Country of Manufacture
	European Conformity marking
	European representative
	Australian sponsor
	Swiss authorized representative

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