## Instructions for Use – EXEL SecureTouch Safety Hypodermic Needles

| REF   | Model                                                | Color Code  |
|-------|------------------------------------------------------|-------------|
| 27403 | EXEL SecureTouch Safety Hypodermic Needle 25G X 5/8" | Orange      |
| 27405 | EXEL SecureTouch Safety Hypodermic Needle 25G X 1"   | Orange      |
| 27406 | EXEL SecureTouch Safety Hypodermic Needle 25G X 1½"  | Orange      |
| 27408 | EXEL SecureTouch Safety Hypodermic Needle 23G X 1"   | Deep Blue   |
| 27416 | EXEL SecureTouch Safety Hypodermic Needle 21G X 1½"  | Deep Green  |
| 27417 | EXEL SecureTouch Safety Hypodermic Needle 20G X 1"   | Yellow      |
| 27420 | EXEL SecureTouch Safety Hypodermic Needle 18G X 1½"  | Pink        |
| 27426 | EXEL SecureTouch Safety Hypodermic Needle 27G X 1½"  | Medium Grey |
| 27427 | EXEL SecureTouch Safety Hypodermic Needle 27G X 1¼"  | Medium Grey |
| 27437 | EXEL SecureTouch Safety Hypodermic Needle 30G X ½"   | Yellow      |

### 1 Device Description

- 1.1 EXEL SecureTouch Safety Hypodermic Needles are single-lumen hypodermic needles with a safety shield mechanism. They consist of a hollow, bevel-edged metal tube (the cannula) that is sharpened at one end and at the other end joined to a female connector (the hub) designed for connection to other devices (e.g., piston syringes, etc.). They are used to inject fluids into, or withdraw fluids from, parts of the body below the surface of the skin. The safety shield mechanism is attached to the hub and is activated after use, covering the needle tip and preventing needlestick injuries and/or re-use of the device. A protective cap is also included to prevent needlestick injuries during transport and handling. These are single-use sterile devices.
- 1.2 Models are differentiated by the needle gauge (the diameter of the sharpened metal tube) and the needle length (the length of the sharpened metal tube). The hub and safety shield are pigmented according to ISO 6009 color-coding to facilitate visible identification of the needle gauge.
- 1.3 Materials:
  - 1.3.1 **Needle Hub / Safety Mechanism:** PP (Ethylene / Polypropylene Copolymer)
  - 1.3.2 **Needle Cannula:** Stainless steel
  - 1.3.3 **Adhesive:** Epoxy Resin
  - 1.3.4 **Lubricant:** Silicone Oil
  - 1.3.5 **Protective Cap:** PP (Ethylene / Polypropylene Copolymer)

### 2 Indications for Use

- 2.1 These devices are intended to be used:
  - 2.1.1 To inject fluids into, or withdraw fluids from, the body.
  - 2.1.2 To administer and/or withdraw (aspirate) fluids/drugs.

### 3 Intended Users / Patient Population

- 3.1 Intended Users: Qualified healthcare professionals in a professional healthcare environment.
- 3.2 Intended Patient Population: The general population, including both adult and pediatric patients of all age groups, as determined to be appropriate by the qualified healthcare professional using the device.

### 4 Clinical Benefit

- 4.1 Safe and reliable delivery of fluids.
- 4.2 Color-coded for easy identification.
- 4.3 Needle remains locked inside activated safety shield, preventing re-use.
- 4.4 Single-use disposables prevent cross-contamination between patients.

### 5 Warnings / Precautions

- 5.1 For single use only. Do not reuse. Reuse of needle can result in contamination or transmit disease. Use once and discard in a designated sharps container in accordance with local safety standards.
- 5.2 Handle with care to avoid needlesticks.
- 5.3 Do not attempt to deactivate the safety device by forcing the needle out of the safety shield.
- 5.4 Do not resterilize. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility, and functional integrity of the device.
- 5.5 Avoid bending or applying excessive force on the needle tube. Excessive force may lead to needle tube breakage. Do not exceed a bending angle of 20°.
- 5.6 If needle is bent or damaged, no attempt should be made to straighten the needle or use the device.
- 5.7 Activation of the safety shield may cause minimum spatter of fluid that may remain on the needle after injection. For best safety, activate the safety shield using a one-handed technique away from self and others.
- 5.8 Not for drug storage.
- 5.9 These devices are not intended for intraocular use. Use of these devices for intraocular administration may cause adverse effects due to the presence of silicone oil.
- 5.10 U.S. Federal Law restricts this device to sale by or on the order of a physician (or other properly licensed practitioner).


### 6 Contraindications


- 6.1 There are no recognized contraindications for the use of these devices.


## 7 Sterility


- 7.1 EXEL SecureTouch Safety Hypodermic Needles are sterilized by EO gas sterilization and delivered sterile in unopened, undamaged packaging.
- 7.2 Single use only. Do not reuse.
- 7.3 Do not resterilize or reprocess.
- 7.4 Inspect each package prior to use. Do not use the device if any seal or cavity is damaged or breached or if the expiration date has been exceeded. Once opened the device must be used or discarded.


## 8 Operating Instructions


- 8.1  Select the appropriate gauge and length of needle for the procedure. Carefully remove the safety needle from the package.


- 8.2  Securely attach the safety needle to a syringe with a push and clockwise twist motion. Move the safety shield away from the needle and toward the syringe barrel. The safety shield will remain in the position you set.

- 8.3  Remove the needle cover carefully to avoid damaging the needle. Perform puncture per standard clinical practice.

- 8.4  After puncture, position shield in preparation for device activation: Using a one-handed technique, push the tab forward with your finger or thumb so that the shield is less than 90 degrees from the needle.

- 8.5  Activate shield: Position the shield approximately 45 degrees to flat surface.

- 8.6  Press down with a firm, quick motion until a distinct audible click is heard.

- 8.7  Visually confirm that the needle is fully engaged under the lock. Dispose of used needles and materials following applicable federal and local safety regulations for sharps disposal.

**Note:** The location of the lock will vary depending on the needle size. Audible click may not be heard on small needle sizes; visual confirmation is required per needle.

## 9 Risks / Side Effects

- 9.1 Use of these devices may expose the user or patient to the following risks and/or side effects:
  - 9.1.1 Misuse of these devices may cause contamination exposing the user or patient to infection, transmission of disease, or bloodborne pathogens. Take care to maintain sterility when using the device and dispose of the device in a sharps container immediately after use.
  - 9.1.2 Misuse of these devices may expose the user or patient to risk of accidental needlestick injuries.
  - 9.1.3 Patients may experience inflammation, hypersensitivity, bruising, bleeding, redness, irritation, or discomfort at the insertion site.
  - 9.1.4 Negative reactions can occur for patients with allergic sensitivity to stainless steel.
  - 9.1.5 Use of these devices for intraocular applications can expose patients to risks due to the presence of silicone oil, including persistent vitreous floaters, inflammation, and increased ocular pressure.

## 10 Storage

- 10.1 Keep away from direct sunlight.
- 10.2 Keep dry.

## NOTICE

© 2026 Exelint International, Co. All rights reserved. EXEL and the EXEL logo are trademarks of EXELINT International, Co.

Any serious incident that occurs must be reported to the manufacturer and competent authority according to local regulations. EXELINT International, Co. disclaims any responsibility for possible consequences resulting from reuse, resterilization or improper use.

**EXELINT International, Co.**  
2500 Santa Fe Avenue,  
Redondo Beach, CA 90278 U.S.A.  
Tel: (310) 649-0707 Fax: (310) 649-1178  
e-mail: [info@exelmed.com](mailto:info@exelmed.com)

FOR U.S.A. INFO CALL: 800-940-3935  
[www.exelint.com](http://www.exelint.com)  
e-mail: [info@exelint.com](mailto:info@exelint.com)

**EMERGO EUROPE**  
Westervoortsedijk 60  
6827 AT Arnhem  
The Netherlands

**EMERGO AUSTRALIA**  
Level 20 Tower II Darling Park  
201 Sussex Street  
Sydney, NSW 2000 Australia

## Graphic Symbol Legend for Medical Device Labeling

|                                                                                     |                                                                                       |
|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|    | Medical Device                                                                        |
|    | Catalog number                                                                        |
|    | Federal law (USA) restricts this device to sale by or on the order of a physician.    |
|    | Consult instructions for use                                                          |
|    | Caution                                                                               |
|    | Sterilized by Ethylene Oxide. Sterility guaranteed if package unopened and undamaged. |
|    | Manufacturer's Lot Number                                                             |
|    | Unique Device Identifier                                                              |
|    | Use by date                                                                           |
|    | Date of Manufacture                                                                   |
|  | For single-patient-use only. Do not reuse.                                            |
|  | Do not resterilize                                                                    |
|  | Device possesses a reuse prevention function.                                         |
|  | Do not use if package is damaged                                                      |
|  | Keep away from sunlight                                                               |
|  | Keep Dry                                                                              |
|  | Non-Pyrogenic                                                                         |
|  | Single Sterile Barrier                                                                |
|  | Protective Packaging Outside of Sterile Barrier System                                |
|  | Manufacturer                                                                          |
|  | Country of Manufacture                                                                |
|  | European Conformity marking                                                           |
|  | European Representative                                                               |
|  | Australian Sponsor                                                                    |

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